



ACC GR 612G-50/616G-60

WHY ACC

- ✓ ACC offers products from the **top-class** activated carbon manufacturers.
- ✓ ACC ensures a **one-stop service** for various demands from customers.
- ✓ **Double quality guarantee** from both manufacturer side and ACC side after production and before delivery.
Approved by many authorized institutions from different fields such as **NSF, BV, HALA, Kosher**.etc.

ACC GR series is coconut-shell-based granular activated carbon designed for gold recovery applications.

It can help to achieve high gold loadings and the relative ease with which they can be stripped enables gold processors to obtain the highest possible yields from ores with high silver/gold ratios, carbonaceous ores, and ores containing other metals that complex with cyanide, for example, copper and nickel.

This product can be effectively re-used in the gold recovery process after being reactivation.

Package: Packed in 500kg jumbo bag. Other packaging is available upon your request.

Standards: All the above test methods are performed using ASTM protocol for activated carbon.

SPECIFICATIONS

GR 612 G-50

| | |
|---------------------|---------------|
| CTC | 50 min/60 min |
| Ash by weight | 4%-5% max |
| Moisture by weight | 5% max |
| Hardness Number | 98% min |
| Platelets by weight | 3% max |
| Particle Size | 6x12 |

GR 616 G-60

| | |
|---------------------|---------|
| CTC | 60 min |
| Ash by weight | 4% max |
| Moisture by weight | 5% max |
| Hardness Number | 98% min |
| Platelets by weight | 3% max |
| Particle Size | 6x16 |