



## ACC GR BA 612/616G

### WHY ACC

- ✓ ACC offers products from the **top-class** activated carbon manufacturers.
- ✓ ACC ensures a **one-stop service** for varies demands from customers.
- ✓ **Double quality guarantee** from both manufacturer side and ACC side after production and before delivery.
- ✓ Approved by many authorized institutions from different fields such as **NSF, BV, HALA, Kosher**.etc.

This series is coconut-shell-based granular activated carbons designed for gold recovery applications.

It can help to achieve high gold loadings and the relative ease with which they can be stripped enables gold processors to obtain the highest possible yields from ores with high silver/gold ratios, carbonaceous ores, and ores containing other metals that complex with cyanide, for example, copper and nickel.

This product can be effectively re-used in the gold recovery process after being reactivation.

**Package:** Packed in 500kg jumbo bag. Other packaging is available upon your request.

**Standards:** All the above test methods are performed using ASTM protocol for activated carbon.

#### SPECIFICATIONS

##### GR BA 612/616G-01

Butane Activity by weight	17.5 min
Ash by weight	4% max
Moisture by weight	5% max
Hardness Number	98% min
Platelets by weight	5% max
Particle Size	6x12 or 6x16

##### GR BA 612/616G-02

Butane Activity by weight	23.3 min
Ash by weight	4% max
Moisture by weight	5% max
Hardness Number	98% min
Platelets by weight	5% max
Particle Size	6x12 or 6x16