



ACC GR K03 for Gold Recovery

WHY ACC

- ✓ ACC offers products from the **top-class** activated carbon manufacturers.
- ✓ ACC ensures a **one-stop service** for varies demands from customers.
- ✓ **Double quality guarantee** from both manufacturer side and ACC side after production and before delivery.
- ✓ Approved by many authorized institutions from different fields such as **NSF, BV, HALA, Kosher**.etc.

This type of coconut shell based activated carbon for gold extraction ,made of original coconut shell, is activated in SLEP activation furnace with optimum technological conditions and the most advanced testing method. The product has the features of highly developed pore structure, large specific surface area, good hardness, excellent adsorption capacity, absorption speed, easy desorption, repeated regeneration etc.

It is mainly used in extracting gold by CIL process.

Package: Packed in 25kg bags. Other packaging is available upon your request.

Standards: All the above test methods are performed using ASTM protocol for activated carbon.

SPECIFICATIONS

GR K03 E

Iodine number (mg/g)	1000 min
Methylene blue (mg/ml)	9 min
Hardness (%)	97 min
Moisture (%)	10 max
Ash (%)	5 max
Particle size (mesh)	90% min pass through 6-12 mesh
PH	7
CTC (%)	55 min
Bulk density (g/ml)	0.40—0.55
Specific surface area (m2/g)	1100 min